

Combining Robust Expectation Maximization and Mean Shift Algorithms for Multiple Sclerosis Brain Segmentation

D. Garcia-Lorenzo^{1 2 3} S. Prima^{1 2 3} D.L. Collins⁴
D.L. Arnold⁴ S.P. Morrissey^{1 2 3 5} C. Barillot^{1 2 3}

1 INRIA, VisAGeS Unit/Project, IRISA, Rennes, France

2 University of Rennes I, CNRS IRISA, Rennes, France

3 INSERM, U746 Unit/Project, IRISA, Rennes, France

4 Montreal Neurological Institute and Hospital, McConnel Brain Imaging Centre, Montreal, Canada

5 Department of Neurology, University Hospital Pontchaillou, Rennes, France

Summary

The logo for VisAGeS features a stylized white dome or canopy structure with a grid of lines, set against a light blue background. Below the dome, the text "VisAGeS" is written in a bold, orange, sans-serif font.

VisAGeS

- Introduction
- Method
- Validation
- Conclusions

Introduction

● Global methods

- Use of the whole image information
- **Expectation-Maximization (EM)**
- Limit: little or none spatial information employed

● Local methods

- Use of local information to create regions
- **Mean Shift (MeS)**
- Challenge: How to merge regions to obtain the final segmentation

- Theory

- Mean Shift (MeS)

- Robust Expectation Maximization (REM)

- REMMeS


Mean Shift

- Non-parametric technique for probability density gradient estimation

$$m_{b,g}(x) = \frac{\sum_{i=1}^n x_i g\left(\left\|\frac{x-x_i}{b}\right\|^2\right)}{\sum_{i=1}^n g\left(\left\|\frac{x-x_i}{b}\right\|^2\right)} - x$$

d = dimension
b = bandwidth
k(x) = kernel profile
g(x) = -k'(x)
n = number of pixels

- A point x will arrive always to a point with $m_{b,g} \approx 0$, called mode. All points arriving to the same mode m will form a region.



Mean Shift for 3D images

- Clustering technique
 - $d = m$ (number of MR sequences)
- Joint spatial-intensity domain
 - $d = m + 3$ (spatial dimensions)
 - Integration of spatial information

$$K_{b_s, b_r}(x) = K\left(\frac{x^s}{b_s}\right) K\left(\frac{x^r}{b_r}\right)$$

Robust Expectation-Maximization

- 3-class Finite Multivariate Gaussian Mixture Model
- Modified Expectation-Maximization algorithm (mEM)
 - Trimmed Likelihood (Neikov et Al. 2006)

$$TL = \sum_{i=1}^{n-h} f(x_{\nu(i)}; \Theta)$$


- Ordering function

$$f(x_{\nu(1)}; \Theta) \geq f(x_{\nu(2)}; \Theta) \geq \dots \geq f(x_{\nu(n)}; \Theta)$$

- In our experiments $h=n/10$


REMMeS

VisAGeS




REMMeS

VisAGeS



REMMeS

VisAGeS



REMMeS

VisAGeS

MR sequences

Normalization

mEM

Normalization of image intensity
to obtain WM variance of 1.0 in
all sequences


Classification

MS Rules

Output
Segmentation


REMMeS

VisAGeS




REMMeS

VisAGeS




REMMeS

VisAGeS




REMMeS

VisAGeS



REMMeS

VisAGeS



Validation : Methods

- Four similar approaches :
 - A. Proposed algorithm
 - B. As A. but replacing mEM with the classical EM
 - C. No MeS regions, each voxel is independently classified with the mEM.
 - D. As C. but replacing mEM with the classical EM

- Synthetic images (T1-w, T2-w and PD-w)
 - MS brain with moderate lesion load from Brainweb
 - 3% noise (n)
 - 0%, 20% and 40% of inhomogeneity (rf)
- Real images (T1-w, T2-w and PD-w)
 - 7 patients
 - Denoised, inhomogeneity correction, normalized in the stereotaxic space, Skull-stripped.
 - Lesions were manually segmented by an expert.


	BW n3rf0	BW n3rf20	BW n3rf40	Average Real
D	0,79	0.80	0.78	-
C	0,72	0.77	0.41	0.52 ± 0.07
B	0,87	0.84	0.79	-
A	0,87	0.85	0.63	0.55 ± 0.05

- Brainweb noise 3% and 0% inhomogeneity:

- Rousseau et al., ISBI'08: 0.63
- Freifeld et al., ISBI' 07: 0.77
- Van-Leemput, TMI'01: 0.80
- Our method: 0.87

Results

- Example of classification with mEM vs. classical EM





VisAGeS


Results

T1-w

T2-w

PD-w

Lesions



D. EM

C. mEM


B. EMMeS

A. REMMeS

Conclusions

- Presented new method combining global and local information
 - Importance of robust Expectation-Maximization
 - MeS will keep the information of the contours
- Future works:
 - Reduce the execution time
 - Automatic adjustment MeS parameters
 - Improve validation (more experts, more methods)

Conclusions



The logo for VisAGeS features a stylized white dome structure with a grid of lines, set against a light blue background. Below the dome is a dark blue silhouette of a landscape with rolling hills.

VisAGeS

Thank you very much for
your attention

Supported by Arsep